



Capsule Filling Machine



Capsule Filling Machine is a versatile lever operated semi automatic machine for various operations with lever order of:

- 1. Plates for putting of capsules by hand.
- 2. Removal of Caps.
- 3. Filling of powder by dispensing of brush pressing.
- 4. Fixing of caps and selling.
- 5. Injection of capsules contact parts of stainless for 100 capsules of size 0.
- 6. Capacity of 200 Capsules on order.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com