



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Spirometer



6 ltr. Capacity with chain compensates counter balance to flat through pulley calibrated to denote volume inlet & outlet tubes stainless steel chamber with corrugated rubber tube and mouth piece,

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com