



Abbe Refractrometer



Abbe refractrometer : for accurate measurement of index of refraction of transparent liquids, solids & solutions. For use in sugar, ghee & oils mills pharmaceutical and chemical laboratories, jam factories and distilleries etc. Measuring range 1.300 to 1.700. Measurement of sugar 0% to 95% complete. designed for accurate and quick examination of the refractive indices & mean dispersion of liquid, solids & oil. nd range : nd 1.300 to nd 1.700 in step of 0.001 accuracy : 0.0002 sugar percentage : 0 to 95 % in steps of 1 division accuracy : + 0.5 % temperature : 0-50 o c.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com