



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Analgesiometer Tail Flick Model



Radiant heart or tail flick type:

Analgesiometer Tail Flick Model provides pain stimulus by heated nichrome wire in a rat's tail to determine analgesic effect of drugs. Complete with one spare element to work on 220 a. C.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com