



Auto Colorimeter



Colorimeter are very useful for analyzing chemicals & blood. Features consist of :
Accurate and direct readability of absorbing levels Rugged, Reliable & easy operations
Consists of a long lasting photo diode detector 1 ml sample size
Battery/Mains operated Perfect durability .

Silent Features:-

- Accurate, Direct Readout of Absorbance
- Rugged, Reliable & Easy to Operate
- Long Life Photo Diode Detector
- 1 ml Sample Size

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Std. glass Filter	8 No.
Minimum Volume	1 ml
Display	3 Digit LED
Range	400-700 nm
Output	OD (0 to 1.99)
Resolution	OD: 0.01
Accuracy	± 0.5% FSI
Detector Photocell	Photodiode
Light Source	3.5V, 100mA
Power Supply	230V, ± 10% AC, 50Hz
Filters	8 Filters- 400, 420, 480, 500, 520, 540, 620, 680 nm
Battery Backup Facility	Optional
Dimension	195X230X95 mm(LXBXH)
Weight	3KG

Accessories:-

- Test Tube (Set of 5)
- Dust Cover
- Instruction Manual

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, Mobile: +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com