



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



## Auto Conductivity-Temp. Meter



- Microprocessor Based
- Highly stable & Accurate
- Auto Ranging Facility
- Cell Constant Measurement and Adjustment Facility
- 16X2 LCD Display
- Automatic Temperature Compensation Facility

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)



## Technical Specification:-

Display	Alphanumeric LCD 16X2
Measurement	Conductivity & Temperature
Range	Five Ranges: - 0-200 $\mu$ S/cm, 0-2 $\mu$ S/cm, 0-20 $\mu$ S/cm, 0-1000 $\mu$ S/cm
Range Selection	Automatic in 5 Auto Ranges
Accuracy Cond. Temp.	$\pm$ 0.5% FS, $\pm$ 1 Digit Temp. $\pm$ 0.1 $^{\circ}$ C $\pm$ 1 Digit
Temperature	Automatic 0-50 $^{\circ}$ C
Cell Constant	Adjustable
Conductivity Cell	Platinum DIP Type
Resolution Cond. Temp	0.1 S/cm, 0.1 $^{\circ}$ C
Power Supply	12V DC Adapter (230V $\pm$ 10% AC, 50 Hz)



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



## Accessories:-

- Conductivity Cell
- Temperature Probe
- Operation Manual
- Dust Cover
- Conductivity Cell Stand
- 12V DC Aadapter

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)