



Automatic Tissue Processor



- All 12 station time can be individually programmed using a pin-up timer
- In-built safety device to protect the tissues in case of power failure
- Imported motor and sturdy gear assembly ensures smooth and silent shifting of tissue baskets
- Thermostatically controlled wax baths with temperature range of ambient to 80 deg C.
- Continuous rotation of tissue baskets ensures better tissue processing

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



TTRIBUTE	M111	M112	M212
Operation	Automat ic	Automatic	Automatic
SS Basket	1	2	2
Basket Rotor	1	2	2
24 hours calibrated discs	6	6	6
Beaker Capacity L	1	1	2
Beakers	10	10	10
Tissue Cassettes Capacity	24	50	100
Tissue Cassette Size (dia)	36mm	36mm	36mm
Safety Device	Yes	Yes	Yes
Cutting Pliers	Yes	Yes	Yes
Wax Baths	2	2	2
Wax Bath Capacity L	1	1	2
Temperature Range deg C	0-85	0-85	0-85
Temperature Cutout deg C	60	60	60

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Vacuum System	Optional	Optional	Optional
Fume Control System	Optional	Optional	Optional

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com