



Binocular Microscope (Led)



- LED LAMP ILLUMINATION
- COAXIAL MECHANICAL STAGE
- BINOCULAR HEAD INCLINED AT 45° & MOVEABLE AT 360°
- CONDENSER ABBE N.A. = 1.2 WITH IRIS DIAPHRAGM & FILTER.
- SEPARATE COARSE & FINE FOCUSING KNOBS
- SPECIAL SOCKET & CORD

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



- OBJECTIVES: 4X , 10X , 40X , 100X (OIL IMMERSION)
- EYEPIECES : 10X WF
- PACKED IN DUST COVER & THERMO PACK

OPTIONAL

- BATTERY BACKUP
- IMPORTED OBJECTIVES

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com