



## Cadaverous Injector



- Cadaverous Injector machine is made of Stainless Steel
- For injecting formaldehyde solution in cadaverous at much higher speed than normal gravity process.
- Unit is fully covered mounted on a portable trolley having four castor wheels for easy movement.
- Unit consists of one air compressor which is connected with a stainless steel tank of 10 liters capacity meant for storing and injecting the solution.
- Incorporated with latest Reciprocatory pressuring unit for efficient performance.
- The machine is Completely Noiseless and Vibration free which justifies the Superiority of our machine from the other brands available.
- Tank is fitted with a safety valve, pressure gauge and rubber tubing having provision for injection.
- Workable on 220 V/50 Hz , Single Phase

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)