



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Clarity Test Equipment for Ampoule



Clarity test equipment for ampoule clarity test equipment for ampoule milky white and black screen fitted on a base board with diffused light arrangement overall dimensions are 23"x15"x18 high. A stainless upright adjustable slide arm with clamp magnifier, main cord and plug.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com