



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Continuous Passive Motion System



continuous-passive-motion-system

Electronic controls enables for range of motion

Adjustable flexion angle from 0 degree to 110 degree

0 to 30 minutes countdown timer

Patient safety switch to stop and reverse the motion

Manual up/down: Manual flexion and release

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Technical Specification of Continuous Passive Motion System

- Mains Input: 20V AC 50Hz
- Power: 20 watt
- Time: 0 to 30 minutes
- Flexion Angle: Adjustable from 0o to 110o
- Delay timer: 0 to 10 second
- Weight:9 kg

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com