



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



## Cook's Pole Climbing Response Apparatus



Built in solid state buzzer and stimulator to provide electrical shocks of 400V in pulsating rates 0.2mA a frequency of 5per sec, for a duration controlled manually or by built in 30 second timer. Output available for recording on kymograph or polygraph complete to work on 220V AC.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)