



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Deep Freezer (Horizontal – 20 C) - Digital



Puff insulated to prevent thermal losses, very low noise plastic body solid door & CFC free heavy duty Compressor, works on 220 V AC / 50 Hz Mains

Temp. Range : (- 20°C)

Total Gross Capacity (L) 100 Ltrs. / 150 Ltrs. / 300 Ltrs.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com