



Digital Dissolved Oxygen Meter



We are frontrunners in offering high quality Digital Dissolved Oxygen Meter. Our R & D team makes sure that our products are regularly tested on internal quality parameters. Certified vendors provide us with raw materials and international quality norms are adhered to during manufacture.

Features:

- - Reliable Instrument
- - Accurate D.O. & Temperature Measurements in the laboratory
- - Available with customized solution
- - Durability and longevity
- - Standard prices
- - Range of design

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Technical Specifications

Display	3½ Digit Bright Red LED
Range	0-20 ppm
Temperature Compensation	Manual 0-50°C
Resolution	0.1 ppm
Accuracy	± 0.1 ppm
Power Supply	230V, ± 10% AC, 50 Hz
Dimension	280X180X88 mm (LXBXH)
Weight	3 KG

Accessories:-

- Gold/Silver Amperometric
- D.O. Probe
- Set of DO Membranes
- Dust Cover
- Operation Manual.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com