



Digital Ph Meter



Features Consist Of High Stability And Accuracy Auto Facility Slope Control 3.5 Digit Led Display Various Available Designs Properly Tested.

Features: -

- High Stability and Accuracy
- Slope Control
- 3.5 Digit LED Display
- Properly Tested
- Original Laboratory Instruments Manufacturer in India

Technical Specification:-

Display	3.5-digit LED
Resolution	0.01 pH, 1 mV

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, Mobile: +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Range	0 to 14.00 pH
Calibration	2 Point (4pH, 7pH)
Accuracy	± 0.01 pH, ± 1 mV
Temperature	Manual 0 to 100 Degree Celsius
Slope Correction	80 to 120%
Power Supply	230V, $\pm 10\%$ AC, 50 Hz
Dimension	280X180X88 mm (LXBXH)
Weight	3 KG

Accessories: -

- Combination pH Electrode
- Buffer Tablets (4pH, 7pH)
- pH Electrode Stand
- Instruction Manual,
- Dust Cover

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com