



## Dissection Table



- Strong Mild Steel tubular Frame Work.
- Top is covered with SS Sheet, with Centre Slope.
- Drainage Hole connected with Waste water outlet Pipe. Provision of SS Bucket and removable Hand rest.
- Mounted on Rubber Stumps. Size: - 73"L x 24"W x 36" H.
- Pre Treated and epoxy Powder coated finishing.
- Also available in half size or in small size.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)