



Distillation Apparatus



The unit is recommended for making pyrogen free distilled water for routine laboratory requirements.

- Made out of Stainless sheet, Heavy Duty.
- A wall- mounting clamp is provided with the unit.
- Supplied with cord and plug to work on 220/230 volts A.C. supply.
- Heaters are of immersion type

Capacity In Ltrs.

- 2 Litre/hrs
- 4 Litre/hrs
- 6 Litre/hrs

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



- 8 Litre/hrs
- 10 Litre/hrs

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com