



## Dry Heat Sterilizer



Dry Heat Sterilizer are designed to suit various applications in growing field of Medical ,Research Laboratories of Hospitals for sterilization. Suitable to work on 220/230V AC supply.

**CONSTRUCTION:-** Double walled construction with inside Stainless Steel Chamber having 60-70mm glass wool insulation & outer wall is made from MS with powder coated. Door is fitted with heavy chrome plated hinges & spring & roller type latch catch. Provided with adjustable steel wire mesh tray for better air circulation.

**TEMPERATURE RANGE:** Ambient to 250°C  $\pm$  1° C. Temperature controlled by MICROPROCESSOR based Digital Temperature Controller or By Hydraulic Thermostat.

**HEATING:** Heating elements are made from high grade (Kanthal ISI) Resistance wire.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.apexambala.com](http://www.apexambala.com)



Heaters are fitted/wound along all sides & bottom. Complete with 2 pilot lamp cord & plug pin etc.

AVAILABLE SIZE		
SIZE INSIDE CHAMBER		
W X D X H	Ltrs.(App.)	No. Of Trays
300 x 300 x 300 mm (12" X 12"X 12")	28	2
355 x 355 x 355 mm (14" X 14"X 14")	45	2
450 x 450 x 450 mm (18" X 18"X 18")	95	2
450 x 450 x 600 mm (18" X 18"X 24")	128	2
600 x 600 x 600 mm (24" X 24"X 24")	227	2
600 x 600 x 900 mm (24" X 24"X 36")	340	2

### **Optional**

- Automatic Holding Timer & Alarm
- Circulation Fan

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, Mobile: +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: [www.apexambala.com](http://www.apexambala.com)



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



*We Can Design Any Size According To Our  
Customer Requirements.*

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.apexambala.com](http://www.apexambala.com)