



Dry Bath Incubator



- Small volume and light weight finish
- Allows for high usage convenience
- 16X2 LCD Display of set and actual time, temperature.
- Comes with built-in temperature calibration function
- automatic fault detection and buzzer alarm function.
- having built-in over-temperature protect in device for safe and stable working.
- Allows easy cleaning and disinfection.

Other details:

- Model comes with PCR thermal cover technology with evaporation prevention function
- Allows simultaneous display of set and actual time, temperature
- Comes with built-in temperature calibration function
- Automatic fault detection and buzzer alarm function
- Having built-in over-temperature protection device for safe and stable working
- Allows easy cleaning and disinfection

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Technical Specifications

Display	16X2 Line Alpha Numeric LCD
Display Format	Menu Driven
Setting	By Soft Touch Keys
Heating Devices	Test Tube Dry Bath
Temp. Sensing Devices	RTD/PT100/Digital Semiconductor
Range	5 Degree Celsius above ambient to 50 Degree Celsius
Resolution	1 Degree Celsius
Accuracy	+/- 1 Degree Celsius
Calibration	Manual (if Required)
Block Size	24 X 12 mm Tube holders
Block Material	Rock Solid Aluminium Block
Power Supply	230V, +/-10% AC, 50Hz

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Accessories:-

- Main Lead
- Dust Cover
- Instruction Manual

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com