



## Fully Automatic Colorimeter



Colorimeter are very useful for analyzing chemicals & blood. Features consist of :  
Accurate and direct readability of absorbing levels Rugged, Reliable & easy operations  
Consists of a long lasting photo diode detector 1 ml sample size  
Battery/Mains operated Perfect durability .

### **Silent Features:-**

- Accurate, Direct Readout of Absorbance.
- Rugged, Reliable & Easy to Operate.
- Long Life Photo Diode Detector.
- 1 ml Sample Size.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Range	400 -700 nm
Accuracy	± 0.01 O.D.
Stability	± 0.02 O.D. Per Hour
Light Source	Highly Accurate Laser Source
Volume	1 ml. Minimum
Dimension	195X215X95 mm (LXBXH)
Battery Backup	6 V DC Re-chargeable Battery (optional)
Display	LED (7 Segment Display)
Display Model	Absorbance & Wavelength
Power Supply	230V, 10% AC, 50Hz
Filters	Built in 9 Digital Filters 50,000 hours Life
Filter Wavelength	400 nm, 430 nm, 480 nm, 500 nm, 520 nm, 540 nm, 580 nm, 620 nm, 680 nm
Weight	1.5 Kg. Approx.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)