



## Hand Operated Tablet Machine



Tablet Machine Hand Operated is a single punch for making tablets upto 1000mg. as per size of the die and punch 6.9 and 12mm. Automatic Compression and ejection device with each Revolution of the wheel drive, automatic feed and removal of tablets. Tablet Machine Hand Operated built on a heavy cast iron body with acentric pressure system stainless steel hopper provided and connected by a rubber tube, easy setting for thickness and hardness of tablets with nuts and lock nuts and set of two spanners. Provided with necessary tool die. Heavy duty.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)