



## High Precrsion Water Bath



- Double walled, Inner chamber made of Stainless steel.
- Outer chamber made of M.S. Duly Powder Coated.
- Gap between the two walls is filled with Glass wool.
- Temperature is controlled by : DIGITAL TEMPERATURE INDICATOR cum CONTROLLER from ambient 5°C to 95°C.
- Stirrer is provided for continuous stirring
- Front glass window on one side for visibility of the sample.
- Workable on 220 Volts 50 Hz. AC mains.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.apexambala.com](http://www.apexambala.com)