



Incubator Shaker



Castor wheel mounted cabinet is double walled with inner chamber of stainless steel and outer casing of M.S. sheet duly painted / S.S. The gap of 75mm between the outer and the inner wall is filled with special grade glass wool to prevent thermal losses. Back of the unit is fitted with air circulation fan for maintaining the temperature uniformly throughout the chamber. Heating elements of nichrome are placed in the air path and are duly insulated from the body.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.apexambala.com



(A) INCUBATOR SHAKER (Orbital Shaker Heating Only)

- With Variable speed and Digital Temperature Controller
- Inner Stainless Steel & Outer M.S. with powder coated.
- Temperature range ambient to 70°C.

Specially designed to hold :

- (i) To Hold 9 flasks of 250 ml OR Mix
(Size of platform 450 x 450 mm)
- (ii) To Hold 12 flasks of 250 ml OR Mix
(Size of platform 600 x 600 mm)

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.apexambala.com