



Lab Centrifuge



These Centrifuges are ideal for growing routine application in biochemical and clinical labs . Unit has digital speed indicator and 0-15 minute digital countdown timer.

Salient Features:

- Stepless speed regulator.
- Safety Lid interlock to prevent cover opening during centrifugation.
- Digital speed meter and 0-15 minutes digital countdown timer.
- Press Switch for quick acceleration to full speed.

Technical details:

Max speed 4000-4500 R.P.M.

With Rotor Head 4 X 50 ml or 8 X 15 ml

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com