



Lab Stirrer



The Motor stirrer having maximum speed upto 3000 RPM unit is clamped to a chrome plated rod fitted on a heavy cast iron base. The height and angle of the stirrer can be adjusted according to the requirement. Variable speed motor with speed regulator is fitted with S.S. stirrer rod to give smooth and noiseless operation. Available with or without Drill Chuck.

Available in following models:

- For 1 & 2 Ltr. (1/20 H.P. Motor)
- For upto 5 Ltr.(1/8 H.P. motor)
- For upto 10 Ltr. (1/4 H.P. Motor)

Optional

- Digital RPM Indicator

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com