



Magnetic Stirrer With Hot Plate



- With S.S. top hotplate.
- Temp of hotplate up to 70°C.
- It has speed regulator for varying speed of stirring.
- Regulator is provided for controlling the RPM of motor
- Electric supply: 220 - 240 V, 50 Hz single phase
- Capacity 1 ltr / 2 Ltrs. / 5 Ltrs. / 10 Ltrs.

Optional

- Energy/RPM Meter/Digital Temperature Controller is provided.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com