



Medical Microscope (Monocular)



- **LIGHT SOURCE :** HEAVY MODEL WITH BUILT IN SOLID STATE TRANSFORMER FOR 6V/20WATT HALOGEN BULB/INCANDESCENT LAMP or LED.
- COAXIAL MECHANICAL STAGE
- MONUCULAR HEAD & MOVEABLE AT 3600
- CONDENSER ABBE N.A. = 1.2 WITH IRIS DIAPHRAGM & FILTER.
- SEPARATE COARSE & FINE FOCUSING KNOBS
- PECIAL SOCKET & CORD

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



- **OBJECTIVES** : 10X , 40X , 100X (OIL IMMERSION)
- **EYEPIECES** : 10X WF
- **HEAD TYPE** : Monocular
- **PACKED IN DUST COVER & THERMO PACK**

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile**: +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com