



Microprocessor Dissolved Oxygen-Temp. Meter



We are frontrunners in offering high quality Digital Dissolved Oxygen Meter. Our R & D team makes sure that our products are regularly tested on internal quality parameters. Certified vendors provide us with raw materials and international quality norms and adhered to during manufacture.

Features: -

- Highly Stable & Accurate
- 16X2 Line Alphanumeric LCD
- 1000 Samples Storage Facility
- PC Interface Facility
- Auto/Manual Calibration
- D.O. & Temperature Measurements
- Soft Touch Membrane Type Keys

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Technical Specifications

Display	16X2 Line Alphanumeric LCD With Backlit
Range	DO: 0-20 ppm, Temp. 0-100°C
Resolution	DO: 0.1 ppm, Temp.: 0.1°C
Accuracy	DO: ± 0.2 ppm, Temp.: ± 0.2 °C
Temp. Comp.	0-50°C (Auto/Manual)
Key Board	8 Keys, Soft Touch Membrane Type
PC Interface	RS 232
Storage	Up to 1000 Samples
Power Supply	230V, $\pm 10\%$ AC, 50 Hz
Dimension	195X225X65 mm (LXBXH)
Weight	2 KG

Accessories:-

- D.O. Probe
- Temperature Probe.
- D.O. Membranes

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



- 12V DC Adapter
- Operation Manual
- Dust Cover
- RS 232
- Software CD

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com