



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



## Multi Chamber Cooling Cabinet



Castor wheel mounted cabinet is double walled with inner chamber of stainless steel 304 GRADE and outer casing of M.S. sheet duly powder coated. The gap of 75mm between the outer and the inner wall is filled with special grade glass wool to prevent thermal losses. Double walled back of the unit is fitted with two air circulation fans for maintaining the temperature uniformly throughout the chamber. The inner chamber is ribbed for placement of shelves at desired heights. Full view inner glass door permits inspection of specimens, without disturbing the same. Lock and key arrangement is provided in the double walled outer door. Unit is fitted with door operated illumination lamp inside the chamber. The temperature range from 0°C to ambient is controlled by Microprocessor Based P.I.D. Digital Temperature Controller cum Indicator. Refrigeration unit is formed by hermetically sealed compressor and the cooling coils are distributed all round the chamber and lie in the air circulation path. Panel consist of separate indicator lamps for mains, heating and incoming voltage. Supplied complete with 2/3 shelves, cord and plug to work on 220/230V A.C. supply.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



## Optional

- Timer
- Fully Automatic Voltage Stabilizer

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)