



Orbital Shaker



Mild Steel Construction Powder Coated, Stainless Steel Platform With Stainless Steel Spring Loaded Lotus Clamps To Hold Flask, Fitted With Heavy Duty Motor Having Variable Speed. Speed Can Be Controlled Through A Continuously Variable Speed Regulator Knob.

SPEED RANGE: 40-200 RPM

Available With the following sizes of platform Size, Capacity And Number Of Flasks.

SIZES : 49 Flasks Of 25ml Or 36 Flasks Of 50ml Or 25 Flasks Of 100ml Or 9 Flasks Of 250ml Or 500ml

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.apexambala.com