



Paraffin Wax Bath



- Double walled outer made of M.S. sheet duly painted / S.S.
- Inner stainless steel sheet.
- The bath consists of 2 pilot lamps & on-off switch.
- Supplied with immersion heaters.
- Temp. Controlled by Thermostat from ambient temp.+5°C to 95°C ± 1°C

Capacity

- 2 Ltr.
- 3 Ltr.
- 4 Ltr..

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com