



Portable Autoclave



This is suitable for general laboratory & field sterilization of instruments & dressing etc. Lid is provided with dial pressure gauge. Double safety valve & steam release valve. Body is made of Aluminum/Powder Coated or Stainless Steel with no joints.

Mode of Heating

Electrically heated by immersion heaters, to work on 220 V/50 Hz. Or By Gas, Burner etc..

SIZE (approx.)	SIZE (approx.)
225 x 225 mm	10 liters
265x290 mm	16 liters
300x300 mm	22 liters
300 x 350 mm	25 liters

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Optional

- Dressing Drum
- Manual Timer
- Digital Timer

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com