



Research Polarimeter



Used in medical, pharmaceutical and chemical research laboratories. Consist of a imported half shade polariser, analyser, 100mm. Diameter circular scale divided in degrees 0-360 deg. X 1 deg. And a sugar scale divided -30 deg. To + 130 deg. X 1 deg. Iss, a divided micrometer drum is given to read 0. 05 deg. (3 minutes) and 0. 1 deg. Iss, the scale is properly illuminated for taking readings, the body is very sturdy and separate sodium lamp is stand is provided. Complete with polarimeter tube, bubble traps at the center, 200mm. Sodium lamp with stand and transformer.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com