



## Rota Rod Apparatus



Rotarod Apparatus engineered by our company comes with digital counter for recording muscle tone, neuro proximity, motor co-ordination in mice and rat, balance, etc. Rota Road Apparatus provided by us consists of compartments that can be two, three or four with 75mm width of each compartment. This Rota rod apparatus possess a rotating road of 25mm diameter with the speed of 5, 10, 15, 20 and 25 rpm. Every compartment of these apparatuses has in built time interval counter. We offer two compartment, three compartment and four compartment Rota Rod apparatus. Today, we are regarded as one of the eminent Digital Rota rod Apparatus Manufacturers in Haryana.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.apexambala.com](http://www.apexambala.com)



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Rota Road Apparatus - Digital Counter for determining neuro-pixicity, muscle tone, balance and motor co-ordination in rats and mice. Consisting of two or three or four compartments of 75 mm width each with a rotating rod of 25 mm diameter having speed of 5,10,15,20 and 25 rpm with time interval counters in each compartment

- Two Compartment
- Three Compartment
- Four Compartment

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.apexambala.com](http://www.apexambala.com)