



## Rotary Spencer Microtome



- Precise sectioning of tissues upto 1 micron thickness
- Front located feed indicator system
- Knife made of Imported fine grain steel alloy
- Convenient specimen holder
- Sturdy rigid construction for vibration free operation
- Protective cover for easy maintenance

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)



# Apex Industrial Electronics

*symbol of quality*



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Section thickness range	1-50 microns
Section thickness settings	1 micron increments
Vertical Specimen Stroke	55mm
Max. specimen size	32mm
Specimen Orientation	360 degrees
Knife Angle	0 to 10 degrees
Knife length (standard)	120mm
Knife holder grip	Single
Specimen Holder	Set of three
Disposable blade holder	Optional
Disposable blade	Optional

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)