



Semi Automatic Microtome



- LCD Display with feather touch keypad
- Cutting/Trimming Mode
- Fast Forward /Backward Function
- Increase/Decrease Section Thickness
- Adjustable Retraction
- Motorised Feed Drive
- Finger Protection Guard
- Standard/Universal Cassette Clamp

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Section Thickness Setting	0.5-100 microns
Setting Values	From 0.5-5 micron in 0.5 micron increments
	From 5-20 micron in 1 micron increments
	From 20-60 micron in 5 micron increments
	From 60-100 microns in 10 micron increments
Trimming Thickness	1-600 microns
Specimen Retraction Range	5-100 microns in 5 micron increments
Object Feed/ Horizontal Feed	28mm
Vertical Stroke	60mm
Specimen Horizontal Feed	28mm
Max. Specimen Size	50 x 50 mm

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com