



Single Beam Uv-Vis Spectrophotometer



1. Standard Quality Scientific Product
2. RS 232C Computer Interface
3. Tight and Nice Shape
4. High Wavelength Precision
5. RS-232C Computer Interface

Technical Specification:-

Optical System	CT, Grating Monochromator
Wavelength Range	195-1000 nm
Wavelength Accuracy	± 2 nm

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Spectrum Bandwidth	2 nm
Wavelength Repeatability	1 nm
Photometric Repeatability	0.2% T
Stability	$\pm 0.004A/h$ @500 nm
Stray Light	$\leq 0.3\%$ @220 nm
Light Source	Deuterium & Tungsten Lamp
Work Pattern	T.A.C.E.
Display Area	0-200% T, -0.3 – 3A
Display Modes	4 LCD
Zero Setting	Automatic
Power Supply	AC 220V ($\pm 10\%$), 50Hz
Weight	12 kg

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com