



Slide Cabinet – Wooden



For keeping 75mm x 25mm glass slides in horizontal position in plastic grooves individually. Cabinet frame of 19mm thick board with melamine lamination. Hinged door with lock is lined with velvet. Drawer frame of ABS Plastic fitted with knob, index card holder. Extra drawer in base with printed index sheets.

- Slide Cabinet, 5 drawers of 200 slides : Total Capacity 1000 slides.
- Slide Cabinet, 10 drawers of 200 slides : Total Capacity 2000 slides.
- Slide Cabinet, 25 drawers of 200 slides : Total Capacity 5000 slides.
- Slide Cabinet, 50 drawers of 200 slides : Total Capacity 10000 slides.
- Slide Cabinet, 100 drawers of 200 slides : Total Capacity 20000

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.apexambala.com