



Stereo Microscope (With Zoom)



Featuring a zoom magnification of 7x to 45x and a long working distance of 100mm, Well suited for biology, medicine, scientific research, modern electronics industry and other industrial uses such as high precision applications. Fully coated prism produces brighter images with higher contrast, and the anti-mold design.

HEAD :

- Binocular or Trinocular heads with 45 Degree inclined tubes.
- Both eye pieces with 1: 5 mm diopter adjustments
- Inter pupillary distance adjustable between 55 mm to 75 mm.

OBJECTIVE:

- Zoom Ratio 6.4: 1, 7x to 45x magnifications.
- Field of view from 33 mm to 4.9 mm.
- Working distance 100 mm.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



STANDS:

- Ergonomically designed pillar or rack and pinion stand with LED illumination and 2 object clamps.
- Alloy metal cast, hardened off-white painted.

ILLUMINATION:

- Transmitted and incident 3W LED illuminations with built-in 100-240V power supply.
- Both illuminators can be used simultaneously and the light intensities can be adjusted separately.

OPTIONAL:

- Auxiliary lens 0.5X, 0.75X, 1.5X & 2.0X.
- USB CAMERAS OR CCD cameras.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com