



Stereo Microscope (Without Zoom)



- Flat Sharp Image erect & unreserved .Wide field eye pieces.
- Para focal from high to low power quick changing paired objectives through 900 rotations.
- Binocular tubes inclined at 450. Adjustable interpupillary distance, Binocular head rotation in 3600.
- Smooth rack & pinion focusing slide way tension controlled large knobs.
- Built - in illumination both for incident light & transmitted light transformer built in base 6V-20W lamps or LED Lamps.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.apexambala.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



- OBJECTIVES : 2x & 4x paired mounted or revolving turret, changing from low to high by simple rotation & no exchange of slides.
- EYEPIECE : Extra wide field 10x
- MAGNIFICATION : 20x to 40x
- EYEPIECES : 10X WF
- PACKED IN DUST COVER & THERMO PACK

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile**: +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.apexambala.com