



Student Microscope



Straight monocular body with triple revolving nose piece having safety stopper for slides. Coarse and fine motion with separate knobs.

OPTICAL COMBINATION

Ach. Objectives : 10x and 45x (Spring Loaded)

Hyg. Eye Pieces : H 10x and H 15x

Illumination : Substage N.A. 1.25 two lens Abbe condenser with iris diaphragm focussable with rack and pinion.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com