



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Student Organ Bath Double Unit



A Perspex organ bath with fine thermostatic temperature control, complete with accessories, bath size 8"x4"x5" with leak proof dovetailed joints, fitted an sturdy metal frame.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com