



Tablet Coating Pan



Tablet Coating pan 12" dia Stainless Steel, this Tablet Coating Pan is worked by oil immersed gear box with F.H.P. motor drive. Air blower is operated independently and air is thermostatically controlled to give temperature 30 to 100 oC. Switches for motor, blower, heater and thermostat are fixed on a panel with neon indicators. Also mains lead with 15 Amp plug provided.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com