



## Tablet Disintegration Test Apparatus Double Basket



The test equipment work on up and down movement of tablet in medium of water kept 37<sub>0</sub>C. Built on a steel cabinet which housed a geared motor for the vertical shaft drive the shaft moves up and down. On the upper end of the shaft fitted a horizontal strip with clamp on either side. These clamps hold glass vessels open at both ends and with stainless steel wire gauge with rubber rings at the lower end, which dips in water. Two independently thermostatically controlled hot plate rating 65watt are fixed on the top of the cabinet upon which beakers of 1000ml size are net filled with water. Temperature is maintained at 37<sub>0</sub>C with and accuracy of 10<sub>0</sub>C. Independent switches with neon indicators for motor hot plate are fitted on a panel. Mains fuse with two Digital Temperature meter and plug provided. Work on 220Volt AC mains.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)