



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



Tablet Disintegration Test Apparatus



Tablet-disintegration with single disc with six removable glass vessels have a stainless steel wire mesh bottom, which dip in water in a beaker of 1000ml kept on the hot plate, temperature of which is maintained 37.0_c with an accuracy of $\pm 0.1_c$.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com