



Tablet Hardness Tester (Pfizer type)



A handy tough and direct dial reading type tablet hardness tester working as pair of pliers generates breaking force in single grip.

This tablet hardness tester offers accurate and reproducible results much faster other conventional hardness testers.

- Accuracy : $\pm 0.2\text{kg}$ Least Count : 0.2Kg
- Safe testing Force : 20kg
- Maximum Tablet Diameter : 25mm

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com