



Tablet Hardness Tester (Monsanto Type)



Tablet Hardness Tester (Monsanto type) is a fairly handy instrument for taking quick readings for the hardness of tablets. This Tablet Hardness Tester (Monsanto type) is made of brass turned parts, heavily chrome plated, the scale is 0 to 20kg and the index mark and division is in 1 Kg reading. Scale adjustable for zero setting as per size of the tablets.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com