



Tablet Making Machine Motorized



This tablet making machine (motorized) useful for making tablets upto 100 mg. Size of punch & die of 9 mm. Automatic compression and ejection device has been provided, actuated by revolution of the wheel drive. In this tablet making machine feed and removal of tablets has also been made automatic. Acentric pressure system has been applied for easy setting of thickness and hardness of tablets. A S.S. hopper has been given.

The unit has been mounted on a mild steel fabrication. Complete with one die and punch set made of stainless steel, easy setting for thickness & hardness of tablets with nuts and lock nuts. We also provide with these machines a set of two spanners.

Specification :

- Motor single phase ½ HP
- Maximum capacity 8-10mm

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com



Apex Industrial Electronics

symbol of quality



An ISO 9001:2015, ISO 13485: 2016

CE CERTIFIED



- Machine speed 60 tablets per minute

Pls Note: We can also offer Heavy duty tablet machine motorized capacity (500 mg) with die & punch 12mm on order basis.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com