



## Tissue Floatation Bath



- Temperature Range ambient +5deg C to 80 deg C.
- Interior of the bath made of black Anodized Aluminum.
- Digital temperature Controller-cum-Indicator
- LED display for Set Value & Process Value
- Tissue floatation bath size: minimum 200 mm X 100 mm

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

**Phone Number :** +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

**Email ID:** apexambala@gmail.com, apexambala@rediffmail.com

**Website:** [www.aie-india.com](http://www.aie-india.com)