



Ultrasonic Cleaner



Ultrasonics (Ultrasound) is the sound which human beings cannot hear. Ultrasound utilizes frequencies of 40 kHz. High frequency electrical energy is converted into ultrasound waves by means of ultrasonic Transducers, which are bonded on the base of S. Cleaning Tank. These high frequency sound waves create in the liquid countless, Microscopic vacuum Bubbles which rapidly expand and collapse.

Salient Features :

- Easy to operate & maintain with one piece S.S. Tank.
- Indigenously manufactured with advanced MOSFET technology with Auto-tuning facility.
- Digital Tuning of Transducers with Generators to avoid any frequency shifted even during demanding applications.
- Compact, rugged and highly durable systems.
- Extensively protected electronic circuits mean longer and safer Operations.
- Available in Sizes: 3, 5, 10, 15, 20, 30 litres or any customized size.

5375, Cross Road No. 1, Ambala Cantt-133001 (Haryana) - India

Phone Number : +91 171-4005686, **Mobile:** +91-9416966084, 9896055375

Email ID: apexambala@gmail.com, apexambala@rediffmail.com

Website: www.aie-india.com